

CODICE PRODOTTO
403



FLACONE DA
100 ML E 360 ML

Disponibile in stock e
private label no MOQ

UNOSOFT HY

EYECARE

COMPOSIZIONE

PoliEsaMetilenBiguanide	0.0001%
EDTA bisodico	0.1%
Sodio ialuronato	0.01%
Eccipienti e Acqua depurata q.b	100%

CARATTERISTICHE DEL PRODOTTO

Soluzione contenente sodio ialuronato e deproteinizzanti è ideale per la disinfezione, conservazione, pulizi ed idratazione di ogni tipo di lente a contatto, sia trasparente che colorata.

L'equilibrata formulazione di UnosoftHY garantisce un'azione antibatterica realmente efficace, notevole comfort e una idratazione di lunga durata.

Conservare le lenti a contatto in sodio ialuronato assicura un'ottima vestibilità delle stesse un miglior comfort nelle prime ore di utilizzo.

PUNTI CHIAVE

Complete solution for the care of all types of soft lenses: disinfects, preserves, cleans, and moisturizes.

Can be sold in a bottle with a lens case included.

100% made in Italy.

PRODUCT CODE
403



**100 ML AND 360 ML
BOTTLE**

**Available in stock and
private label no MOQ**

UNOSOFT HY

EYECARE

COMPOSITION

Polyethylene biguanide	0.0001%
Disodium EDTA	0.1%
Sodium hyaluronate	0.01%
Excipients and purified water q.s.	100%

PRODUCT FEATURES

This solution containing sodium hyaluronate and deproteinizing agents is ideal for disinfecting, storing, cleaning, and moisturizing all types of contact lenses, both clear and colored.

The balanced formulation of UnosoftHY guarantees truly effective antibacterial action, remarkable comfort, and long-lasting hydration.

Storing contact lenses in sodium hyaluronate ensures excellent fit and improved comfort during the first hours of use.

KEY POINTS

Complete solution for the care of all types of soft lenses: disinfects, preserves, cleans, and moisturizes.

Can be sold in a bottle with a lens case included.

100% made in Italy.